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References: 1 Product Manual - Fluid management under compression for seven days, PD-741114, 2024 2 Product Manual - Fluid Handling Capacity, PD-740029, 2024 5 Lev-Tov H. A prospective, open, multi-center, interventional, noncomparative clinical investigation to follow the progress of exuding chronic wounds using Mepilex Up, Clinical Investigation Report - Diplo01, Mölnlycke Health Care, PD-737729 Rev 1, 2024. 6 Product Manual - Absorption capacity using inclined plane, PD- 741137, 2024 7 Van Overschelde, P. et al. A randomised controlled trial comparing two wound dressings used after elective hip and knee arthroplasty. Poster presentation at 5th Congress of the WUWHS, Florence, Italy, 2016. 8 Bredow J. et al. Evaluation of Absorbent Versus Conventional Wound Dressing. A Randomized Controlled Study in Orthopedic Surgery. Deutsche Arzteblatt International, 2018. 9 Silverstein P. et al. An open, parallel, randomized, comparative, multicenter study to evaluate the cost-effectiveness, performance, tolerance, and safety of a silver-containing soft silicone foam. Journal of Burn Care and Research, 2011. 10 Gee Kee E.L. et al. Randomized controlled trial of three burns dressings for partial thickness burns in children. Burns, 2014. 11 David F. et al. A randomised, controlled, non-inferiority trial comparing the performance of a soft silicone-coated wound contact layer (Mepitel One) with a lipidocolloid wound contact layer (UrgoTul) in the treatment of acute wounds. International Wound Journal, 2017. 12 Gotschall C.S. et al. Prospective, randomized study of the efficacy of Mepitel on children with partial-thickness scalds. Journal of Burn Care & Rehabilitation, 1998.

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© 2024 Mölnlycke Health Care. Alle Rechte vorbehalten. DWC0204